



WARM REMINDER

Dear Client:

Thank you for choosing COXO Oral Scanner. In order to make the equipment work normally, please pay attention to the following key points:

1. Prior to install the scanning software, please ensure that your computer fulfills the configuration requirements. If the computing power is not enough, it will lead to lagging of scanning and long processing time! The configuration requirements in detail are as follows:

	PC configuration requirement	Note
CPU	Intel Core i7/i9 (10th/11th generation), with base frequency of 3.0GHz or above ; AMD R7 5800H or above, with base frequency of 3.0GHz or above;	Not compatible to Apple CPU
GPU	NVIDIA GeForce RTX 4060 or above	Not compatible to AMD GPU
RAM	16GB RAM or above	
Monitor	Monitor with HD resolution of 1920X1080 or above	
Operating system	Windows 7/10/11	Not compatible to Apple IMAC
USB port	USB 3.0/3.2/3.2 (at least 2)	

2. In order to guarantee the computing power of the computer, please make sure to plug in the power cord and set it to the best performance mode so as to ensure the smooth running of the software.

3. During usage, please ensure that the computer remains connected to WIFI and the signal is stable; otherwise, the software will not function properly.